Nom: _____ Classe: _____ Date : _____

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Partie 1: Lettre 'e'

1. Cut out a small letter "e" from the newspaper. Place it on a slide, put a drop of water on it and cover it with a cover slip. Use a small ruler to measure the size distance across the letter e. Record your observations.
2. Place the slide with the specimen on the stage. Use stage clips to hold it in place. **MAKE SURE THE LETTER "e" IS POSITIONED UPRIGHT UNDER THE MICROSCOPE.**


3. Turn the light on and maneuver the slide into position so that the specimen sits on top of the hole in the stage.
4. With the low power (10X) objective lens, use the course and fine adjusting knobs to view the letter "e" clearly. Create a **scale drawing** of what you see.
5. Move the microscope slide from left to right and note what direction the letter e moves in your field of view.
6. Move the microscope slide up and down and note what direction the letter e moves in your field of view.
7. Take a clear plastic ruler and place it over top of the microscope slide to measure the diameter of your field of view. Record your observations.
8. Estimate the size of the letter e according to how many would fit across the field of view. Record your observations.
9. Move the revolving nose piece to the medium power lens and focus using **ONLY WITH THE FINE ADJUSTMENT KNOB**. Note what happens when you change the lenses.
10. Measure and record the distance across the field of view using the clear plastic ruler.

Partie 2: Montage humide


1. Using the tweezers, peel a thin transparent layer of soft skin from the outside of an onion and cut it into a 0.5 cm x 0.5 cm square.
2. Carefully drop the onion piece onto a clean microscope slide without trapping air bubbles underneath it.
3. Add 1-2 drops of iodine to the onion skin.
4. Place a cover slip on the onion piece by holding one edge of the slip just inside the liquid and holding the slip at a 45 degree angle. Slowly lower the slip into place to push out any air bubbles.
5. View under the microscope and record your observations.

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